

Datasheet

IQCB1 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009657-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant IQCB1.

Immunogen: IQCB1 (NP_001018865.1, 1 a.a. ~ 465
a.a) full-length recombinant protein with GST tag. MW of
the GST tag alone is 26 KDa.

Sequence:

MKPTGTDPRILSIAAEVAKSPEQNVPVILLKLKEIINITPL
GSSELKKIKQDIYCYDLIQYCLLVLSQDYSRIQGGWTTI
SQLTQILSHCCVGLPGEDAEEFYNEILLPSAAENFLVL
GRQLQTCFINAAKAEKDELLHFFQIVTDSLFWLLGGH
VELIQNVLQSDHFLHLLQADNVQIGSAVMMMLQNILQI
NRSKRSMLEINRQKEEEDLKLQLQLQRQRAMRLSR
ELQLSMLEIVHPGQVEKHYREMEEKSALI IQKHWRGY
RERKNFHQQRQSLIEYKAAVTLQRAALKFLAKYRKKKK
LFAPWRGLQELTDARRVELKKRVDDYVRRHLGSPMS
DVVSRELHAQAQERLQHYFMGRALEERAQQHREALIA
QISTNVEQLMKAPSLKEAEGKEPELFLSRSRPVAAKAK
QAHLTTLKHIQAPWWKKLGEESGDEIDVPKDELSIELE
NLFIGGTKPP

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 9657

Gene Symbol: IQCB1

Gene Alias: NPHP5, PIQ, SLSN5