

Datasheet

AMMECR1 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009949-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant AMMECR1.

Immunogen: AMMECR1 (NP_001020751.1, 1 a.a. ~ 296 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MAAGCCGVKKQKLSSSPSGSGGGGASSSSHCSG
ESQCRAGELGLGGAGTRLNGLGGLTGGGSGSGCTLS
PPQGC GG GGGIALSPPPSCGVGTLLSTPAAATSSSP
SSSSAASSSSPGSRKMVSAEMCCFCFDVLYCHLYG
YQQPRTPRFTNEPYALKDSRFPPMTRDELPRLFCSVS
LLTNFEDVCDYLDWEVGVHGIRIEFINEKGSKRTATYL
PEVAKEQGW DHIQTIDSLLRKGGYKAPITNEFRKTIKLT
RYRSEKMTLSYAEYLAHRQH H FQNGIGHPLPPYNHYS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 9949

Gene Symbol: AMMECR1

Gene Alias: AMMERC1