



ZPACK13951035 PLASTIC BOTTLE LIQ 5 L DIN 45 UN X



SPECIFICATIONS

UN-certification :	3H1/X1.5/Y1.9/Z1.9/250/...
Material :	HDPE
Neck Type :	DIN 45
Filling Opening :	45,00 mm
Height :	332,40 mm
Length :	168,00 mm
Wide :	168,00 mm

AVAILABLE ACCESSORIES

Adaptator for dispenser S40 to GL 45 :	217-9271
Key for opening :	84544.001

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343