



BBP™33 Label Printer - Lab and CLP Kit

The way label making should be



The BBP33 Label Printer is a completely hassle free, easy to use PC-based printer that can print signs and labels up to 101.6 mm wide. Smart chip supplies automatically set up the label in your Brady software, eliminating manual calibration and label waste. Drop-in ribbon and label cartridges ensure fast and effortless material changes. The BBP33 Label Printer can print on a full range of laboratory labels and provides you clear and durable labelling to eliminate errors in sample identification and ensure traceability.

It is really that simple and hassle free – exactly the way label making should be.

Solves all your labelling challenges

- Automatic label setup with no calibration
- Fast and effortless material changeovers
- Prints correctly on the first label!

Impressive material lineup for your Laboratory ID

Depending on the customer's needs, the BBP33 Label Printer can be used with a wide range of materials specially designed for Laboratory Identification:

- Vial & Tube
- Vial Top Labels
- Well Plate Labels
- Centrifuge, Eppendorf & PCR Tube Labels
- Conical / Large Tube Labels
- Slide Labels
- Straw Labelling
- Tissue Cassette Labels
- Bottles, Flasks, Dishes & General ID



VWR Lab kit BBP33 UK: Catalogue number 116-0329	VWR Lab kit BBP33 EU: Catalogue number 116-0330
BBP33 printer & software LabelMark	
Bottles & flasks labels B33-617-488	
Vial & tube labels B33-133-461	
Ink ribbon B30-R4300	

BBP™33 Label Printer

Fast supply changeover and waste-free setup, never-before-seen in a PC-based printer.



- Thermal Transfer Printer
- Print Resolution: 300 dpi
- Print speed: 101 mm/sec. (Max.)
- Prints die cut labels from 6.35 mm to 101 mm wide
- Continuous labels from 12.5 mm up to 3 m long
- Built-in bar code symbologies
- Dimensions: 230.6 x 241.6 x 304.8 mm (W x H x D)
- Use with LabelMark™, CodeSoft™ and MarkWare™ Labelling Software or Windows® based driver for 3rd party software use



Automatic label setup with no calibration

The printer “reads” critical information from your cartridges and automatically sets up the correct label format in the printer or your Brady software.

Fast and effortless material changeovers

Drop in the cartridge, lock it and start typing – that’s all you need to do! (anyone can do it!)



Prints correctly on the first label!

No more wasted labels - save labels and money



VWR CLP kit BBP33 UK: Catalogue number 116-0331	VWR CLP kit BBP33 EU: Catalogue number 116-0332
BBP33 printer & software MarkWare	
CLP labels B30-260-7569-CLP1	
CLP labels B30-260-7569-CLP2	
CLP labels B30-261-7569-CLP2	
Continuous white roll B30C-2250-7569-WT	
Ink ribbon B30-R10000	



www.vwr.com