

Datasheet

IFI44 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00010561-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human IFI44 protein.

Immunogen: IFI44 (NP_006408.2, 1 a.a. ~ 444 a.a) full-length human protein.

Sequence:

MAVTTRLTRLHEKILQNHFGGKRLSLLYKGSVHGFRN
GVLLDRCCNQGP TLTVIYSEDHII GAYAEESYQEGKYA
SIILFALQDTKISEWKLGLCTPETLFCCDVTKYNSPTNF
QIDGRNRKVIMDLKT MENLGLAQNCTISIQDYEVFRCE
DSLDERKIKGVIELRKSLLSALRTYEPYGSLVQQIRILL
GPIGAGKSSFFNSVR SVFQGHVTHQALVGTNTTGISE
KYRTYSIRDGKDGKYL PFILCDSLGLSEKEGGLCRDDI
FYILNGNIRDYQFNPMESIKLNHHDYIDSPSLKDRIHC
VAFVFDASSIQYFSSQMIVKIKRIRRELVNAGVVHVALL
THVDSMDLITKGD LIEIERCEPVR SKLEEVQRKLG FALS
DISVVSNY SSEWELDPVKDVLILSALRRMLWAADD FLE
DLPFEQIGNLREEIINCAQGKK

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10561

Gene Symbol: IFI44

Gene Alias: MTAP44, p44